

## Socialization Rack 6.25 sq. ft.

## Catalog # 38399

The Socialization Rack system is constructed of stainless steel.

The rack is pairable both vertically and horizontally, and has two perches. Socialization panels are provided on each side for 2 phases of interaction. The cage door is on hinges, and can be fully opened, it contains an integrated horizontal sliding animal door. Metabolic pans can be used on both cages. There is a push handle at the rear of the rack. A removable food box, and foraging tray are provided.

The back wire panel of the cages allow more ambient light into the cage.

A single handed squeeze pull also acts as a perch. The squeeze is operated by pulling the perch, it glides smoothly and quietly on Tivar Hot Plastic.

There is an access door on the front of the cage for biscuit retrieval or specimen collection.

Waterlines are mounted on the side of the rack for easy access to test the lixits and drains.

All pairing components are completely self contained. The waste pan, upper floor, and lower ceiling telescopes to the rear for vertical pairing.

## **Cage Specifications:**

Each cage is 6.25 Ft.2 Latching Squeeze System Suitable for 250F Autoclave Stainless Steel Water System 5" Colson Casters





