

# Video Socialization Rack

## Catalog # 37746

The Socialization Rack system is constructed of stainless steel. The rack is pairable both vertically and horizontally, and has two 3-bar perches. The cage door is on hinges, and can be fully opened, or removed, it contains an integrated horizontal sliding animal door. Metabolic pans can be used on both cages. There is a push handle at the rear of the rack. A stationary food box, biscuit retrieval door and foraging tray are included.

The divider panels on the Video Socialization Rack are fully removable. You simply pull the panel off of the track and insert the lexan panel for unobstructed observation of the animal. The vertical pairing feature is fully functional with either wall panel. Wall panels are latched in for added security. The rack can pair with other similar units such as our Socialization Rack when the divider panel is in place.

The vertical pairing feature is self contained. No parts are removed, and no storage is needed. The top cage floor and waste pan telescope to the rear of the cage. This allows the waste in the top waste pan to flow to the rear trough in either the open or closed position. Latches secure the floor in place. A mid-section door on the front of the cage allows added security, and is used to lock a metabolic pan in place (optional).

### **Cage Specifications:**

- Each cage is 6.25 Ft.2
- Latching Squeeze System
- Suitable for 250F Autoclave
- Stainless Steel Water System
- 5" Colson Casters

