

Socialization Rack 6.25 sq. ft.

Catalog # 39082

The Socialization Rack system is constructed of stainless steel.

The rack is pairable both vertically and horizontally, and has one perch in each cage. Socialization panels are provided on each side for 2 phases of interaction. The cage door is on hinges, and can be fully opened, it contains an integrated horizontal sliding animal door. Metabolic pans can be used on both cages. There is a push handle at the rear of the rack. A removable food box, and foraging tray are provided.

A single handed squeeze pull also acts as a perch. The squeeze is operated by pulling the perch, it glides smoothly and quietly on Tivar Hot Plastic.

There is an access door on the front of the cage for biscuit retrieval or specimen collection.

Rack has a dry waste systems. The mid section floor and waste pan are removable. The rear half of the floor remains when the waste pan / floor cartridge is removed.



Each cage is 6.25 Ft.2 Latching Squeeze System Suitable for 250F Autoclave Stainless Steel Water System 5" Colson Casters





