

Socialization Rack

6.25 Sq. Ft.

Catalog # 39963

The Socialization Rack system is constructed of stainless steel. This rack is identical to our Socialization Rack #35465. The one differentiating feature being our Removable Telescoping Floor Cartridge.

The rack is pairable both vertically and horizontally, and has one 3-bar perch and one 4-bar perch. Socialization panels are provided on each side for 3 phases of interaction. The cage door is on hinges, and can be fully opened, or removed, it contains an integrated horizontal sliding animal door. Metabolic pans can be used on both cages. There is a push handle at the rear of the rack. A stationary food box, biscuit retrieval door and foraging tray are included.

The Removable Telescoping Floor Cartridge makes the Socialization Rack more versatile. By removing the mid-section floor it allows the Rack to be a Play Cage. A climbing pole can be added for animal enrichment.

The Floor Cartridge allows for the top cage floor and waste pan to telescope to the rear of the cage. This allows the waste in the top waste pan to flow to the rear trough in either the open or closed position. Latches secure the floor in place. A mid-section door on the front of the cage allows added security, and is used to lock a metabolic pan in place (optional).

Cage Specifications:
Each cage is 6.25 Ft.²
Latching Squeeze System
Suitable for 250F Autoclave
Stainless Steel Water System
5" Colson Casters

